

Clinical Evaluation Report

1. Purpose:

In order to verify the clinical performance of the registered test, this clinical evaluation is conducted.

2. Product information:

COVID-19 Antigen Rapid Test Kit(Swab) was produced by Safecare Biotech(Hangzhou) Co.,Ltd., Lot number is COV20081001, valid until August,2022.

3. Sample requirement:

Fresh samples were collected from CDC and validated by PCR.

4. Supporting equipment:

PCR tests are performed on ABI7500.

The test-strips are manually operated and visually interpreted.

5. Clinical evaluation:

The test was preformed by Wei Lihua/Zhen CaiWen at R&D lab.

6. Statistical methods:

		Referencing reagent Test		Т-4-1
		Positive	Negative	Total
Research Reagent	Positive	А	В	A+B
	Negative	С	D	C+D
Total		A+C	B+D	A+B+C+D

Percent Positive Agreement=A/(A+C)*100%

Negative Percent Agreement=D/(B+D)*100%

7. Evaluation indicators:

The total PPA should be no less than 80%.

The total NPA should be no less than 90%.

- 8. The test data: Refer to the Data Sheet.
- 9. Statistical results of the clinical evaluation

		Referencing Method (RT-PCR)		T-4-1
		Positive	Negative	Total
Test-strip	Positive	30	1	31
	Negative	3	52	55
Total		33	53	86

Statistical results

Project	Value	Percentage (95% confidence interval)
Relative Sensitivity (%)	30/33	90.91% (75.67% ~ 98.08%)
Relative Specificity (%)	52/53	98.11% (89.93% ~ 99.95%)
Positive expectation Rate (%)	30/31	96.77% (83.30% ~ 99.92%)
Negative expected Rate (%)	52/55	94.55% (84.88% ~ 98.86%)
Overall Agreement (%)	82/86	95.35% (88.52% ~ 98.72%)



3) Kappa consistency test

According to the literature [5.1], Calculate the Kappa value and standard error; test hypothesis is established for Kappa: H0: k = 0, Kappa value comes from 0 population, H1: k > 0, Kappa value comes from non-0 population, $\alpha = 0.05$.

Project	Value
Kappa Value	0.9005, Good consistency.
Standard Error Se(K)	0.0485
95% Confidence Interval	0.8055 ~ 0.9956
Standard Error SeO(K)	0.108
Test Value Z	Z=8.3614, Probability value P=0.0000
Test Result	P<0.05,refuse H0, Kappa values come from
	populations other than 0.

4) Conclusion

A side-by-side comparison was conducted using the research reagent and referencing reagent. Compare with RT-PCR:

The Relative Sensitivity is 90.91%, the Relative Specificity is 98.111%, the Overall Agreement is 95.35%.

In summary, The study showed that there is a high coincidence rate between the test-strip and RT-PCR, and have the equivalence on the clinical usage.