

MediSafe® Solo

Lite & Simple

MEDISAFE® SOLO SAFETY LANCETS

Sterile, single-use medical devices
intended for capillary blood sampling



Product intended use

MediSafe® Solo safety lancets are sterile, single-use medical devices intended for capillary blood sampling. Intended users are healthcare professionals.

Product Features

Lite


Simple

Lite product design and geometry of the needle provides adequate blood sample

Flat push-button activation with extra ribs for comfortable finger placement

Ergonomic shape with extra side wings for a firm grip

Long protective cap guarantees needle sterility and gives a secure grip for easy twist off removal



Sharps injuries expose healthcare workers to blood-borne pathogens. The Centers for Disease Control and Prevention (CDC) estimates that about 385,000 sharps-related injuries occur annually among healthcare workers in hospitals.

It has been estimated that about half or more of sharps injuries go unreported. Most reported sharps injuries involve nurses, but laboratory workers, physicians, housekeepers, medical students and other healthcare workers are also exposed to injuries.

Sharp injuries can occur before, during or even after use of a sharp and continue to be a serious health and safety hazard.

Sharps injury - facts

“Over 80% of needlestick injuries can be prevented with the use of safer needle devices, which, in conjunction with worker education and work practice controls, can reduce injuries by over 90%!”


HTL-STREFA is committed to manufacturing safety products designed to reduce the risk of exposure to blood-borne pathogens.

Safety lancets, including **MediSafe® Solo**, are an integral component to the sharps injury prevention programs in hospitals, clinics, laboratories, doctor's offices and wherever both patients and professionals need to feel safe.


By using safety lancets needlestick injuries and infections can be minimized.

**The WHO estimates based on CDC Report*


How to use MediSafe® Solo:



Twist off the protective cap and pull it straight out.



Press the **MediSafe® Solo** body firmly against the puncture site and push the button to activate the device.



Gently apply intermittent pressure near the puncture site to obtain the required blood volume. Discard the used lancet into a sharps container in accordance with facility guidelines and local regulations.

MediSafe® Solo

Lite & Simple

MediSafe® Solo is sterile, single-use, push button activated safety lancet intended for capillary blood sampling.

Available in two versions:

Gauge

Penetration depth



MediSafe® Solo

29G needle

1.5 mm



MediSafe® Solo

23G needle

2.0 mm

References:

Directions for safety lancets development based on nurses' needs; Poland, June 2016, HTL-STREFA S.A. data on file
<https://www.osha.gov/needlesticks/needlefaq.html>

Mary Foley, MS, RN, Annemarie T. Leyden, EdD, RN; American Nurses Association – Independent Study Module
Needlestick Safety and Prevention; http://www.who.int/occupational_health/activities/1anaism.pdf

Poland



HTL-STREFA S.A.
ul. Adamówek 7
95-035 Ozorków
info@htl-strefa.pl
www.htl-strefa.com



Find us on Facebook
www.facebook.com/HTLSTREFA

USA

HTL-STREFA, Inc.
3005 Chastain Meadows Pkwy
Suite 300
Marietta, GA 30066

+1 770 528 0410
info@htl-strefa.com
www.htl-strefa.com



520/ALL/3/2018
Made in Poland

CE 0344

CE mark refers to the product
not to the brochure