









Buprenorphine	5	Norbuprenorphine	25
Buprenorphine 3-D-Glucuronide	25	Norbuprenorphine 3-D-Glucuronide	50
<b>COCAINE (COC 1,500)</b>			
Benzoylcegonine	1,500	Cocaethylene	100,000
Cocaine HCl	1200	Ecgonine	150,000
<b>COCAINE (COC 300)</b>			
Benzoylcegonine	300	Cocaethylene	20,000
Cocaine HCl	200	Ecgonine	30,000
<b>COCAINE (COC 200)</b>			
Benzoylcegonine	200	Cocaethylene	13,500
Cocaine HCl	135	Ecgonine	20,000
<b>COCAINE (COC 150)</b>			
Benzoylcegonine	150	Cocaethylene	1,000
Cocaine HCl	120	Ecgonine	15,000
<b>COCAINE (COC 100)</b>			
Benzoylcegonine	100	Cocaethylene	7,000
Cocaine HCl	80	Ecgonine	10,000
<b>MARIJUANA (THC 300)</b>			
Cannabinol	200,000	$\Delta^9$ -THC	100,000
11-nor- $\Delta^9$ -THC-9 COOH	200	$\Delta^9$ -THC	100,000
11-nor- $\Delta^9$ -THC-9 COOH	300		
<b>MARIJUANA (THC 200)</b>			
Cannabinol	140,000	$\Delta^9$ -THC	68,000
11-nor- $\Delta^9$ -THC-9 COOH	120	$\Delta^9$ -THC	68,000
11-nor- $\Delta^9$ -THC-9 COOH	200		
<b>MARIJUANA (THC 150)</b>			
Cannabinol	100,000	$\Delta^9$ -THC	50,000
11-nor- $\Delta^9$ -THC-9 COOH	100	$\Delta^9$ -THC	50,000
11-nor- $\Delta^9$ -THC-9 COOH	150		
<b>MARIJUANA (THC 50)</b>			
Cannabinol	35,000	$\Delta^9$ -THC	17,000
11-nor- $\Delta^9$ -THC-9 COOH	30	$\Delta^9$ -THC	17,000
11-nor- $\Delta^9$ -THC-9 COOH	50		
<b>MARIJUANA (THC 30)</b>			
Cannabinol	20,000	$\Delta^9$ -THC	10,000
11-nor- $\Delta^9$ -THC-9 COOH	20	$\Delta^9$ -THC	10,000
11-nor- $\Delta^9$ -THC-9 COOH	30		
<b>MARIJUANA (THC 25)</b>			
Cannabinol	17,500	$\Delta^9$ -THC	8,500
11-nor- $\Delta^9$ -THC-9 COOH	15	$\Delta^9$ -THC	8,500
11-nor- $\Delta^9$ -THC-9 COOH	25		
<b>MARIJUANA (THC 20)</b>			
Cannabinol	14,000	$\Delta^9$ -THC	6,800
11-nor- $\Delta^9$ -THC-9 COOH	12	$\Delta^9$ -THC	6,800
11-nor- $\Delta^9$ -THC-9 COOH	20		
<b>METHADONE (MTD 300)</b>			
Methadone	300	Doxylamine	100,000
<b>METHADONE (MTD 200)</b>			
Methadone	200	Doxylamine	65,000
<b>METHAMPHETAMINE (MET 1, 000)</b>			
p-Hydroxymethamphetamine	25,000	(±)-3,4-Methylenedioxy-methamphetamine	12,500
D-Methamphetamine	1,000		
L-Methamphetamine	20,000	Mephentermine	50,000
<b>METHAMPHETAMINE (MET 500)</b>			
p-Hydroxymethamphetamine	12,500	(±)-3,4-Methylenedioxy-methamphetamine	6,250
D-Methamphetamine	500		
L-Methamphetamine	10,000	Mephentermine	25,000
<b>METHAMPHETAMINE (MET 300)</b>			
p-Hydroxymethamphetamine	7,500	(±)-3,4-Methylenedioxy-methamphetamine	3,750
D-Methamphetamine	300		
L-Methamphetamine	6,000	Mephentermine	15,000

<b>METHAMPHETAMINE (MET 200)</b>			
p-Hydroxymethamphetamine	5,000	(±)-3,4-Methylenedioxy-methamphetamine	2,500
D-Methamphetamine	200		
L-Methamphetamine	4,000	Mephentermine	10,000
<b>METHYLENEDIOXYMETHAMPHETAMINE (MDMA 1, 000) Ecstasy</b>			
(±) 3,4-Methylenedioxy-methamphetamine HCl	1,000	3,4-Methylenedioxyethyl-amphetamine	600
(±) 3,4-Methylenedioxyampheta mine HCl	6,000		
<b>METHYLENEDIOXYMETHAMPHETAMINE (MDMA 500) Ecstasy</b>			
(±) 3,4-Methylenedioxy-methamphetamine HCl	500	3,4-Methylenedioxyethyl-amphetamine	300
(±) 3,4-Methylenedioxyampheta mine HCl	3,000		
<b>METHYLENEDIOXYMETHAMPHETAMINE (MDMA 300) Ecstasy</b>			
(±) 3,4-Methylenedioxy-methamphetamine HCl	300	3,4-Methylenedioxyethyl-amphetamine	180
(±) 3,4-Methylenedioxyampheta mine HCl	1,800		
<b>MORPHINE (MOP/OPI 300)</b>			
Codeine	200	Norcodeine	6,000
Levorphanol	1,500	Normorphine	50,000
Morphine-3-β-D-Glucuronide	800	Oxycodone	30,000
Ethylmorphine	6,000	Oxymorphone	50,000
Hydrocodone	50,000	Procaine	15,000
Hydromorphone	3,000	Thebaine	6,000
β-Monoacetylmorphine	300	Morphine	300
<b>MORPHINE (MOP/OPI 200)</b>			
Codeine	160	Norcodeine	4,000
Levorphanol	1,000	Normorphine	40,000
Morphine-3-β-D-Glucuronide	800	Oxycodone	20,000
Ethylmorphine	4,000	Oxymorphone	40,000
Hydrocodone	40,000	Procaine	10,000
Hydromorphone	2,000	Thebaine	4,000
β-Monoacetylmorphine	200	Morphine	200
<b>MORPHINE (MOP/OPI 100)</b>			
Codeine	80	Norcodeine	2,000
Levorphanol	500	Normorphine	20,000
Morphine-3-β-D-Glucuronide	300	Oxycodone	10,000
Ethylmorphine	2,000	Oxymorphone	20,000
Hydrocodone	20,000	Procaine	5,000
Hydromorphone	1,000	Thebaine	2,000
β-Monoacetylmorphine	200	Morphine	100
<b>METHAQUALONE (MQL 300)</b>			
Methaqualone	300		
<b>MORPHINE/OPIATE (OPI 2,000)</b>			
Codeine	2,000	Morphine	2,000
Ethylmorphine	3,000	Norcodeine	25,000
Hydrocodone	50,000	Normorphine	50,000
Hydromorphone	15,000	Oxycodone	25,000
Levorphanol	25,000	Oxymorphone	25,000
β-Monoacetylmorphine	3,000	Procaine	50,000
Morphine 3-β-D-glucuronide	2,000	Thebaine	25,000
<b>MORPHINE/OPIATE (OPI 1,000)</b>			
Codeine	1,000	Morphine	1,000
Ethylmorphine	1,500	Norcodeine	12,500
Hydrocodone	25,000	Normorphine	25,000
Hydromorphone	7,500	Oxycodone	12,500
Levorphanol	12,500	Oxymorphone	12,500
β-Monoacetylmorphine	1,500	Procaine	25,000
Morphine 3-β-D-glucuronide	1,000	Thebaine	12,500
<b>MEPERIDINE (MPRD 100)</b>			
Norperidine	100	Meperidine	100

<b>PHENCYCLIDINE (PCP 50)</b>			
Phencyclidine	50	4-Hydroxyphencyclidine	25,000
<b>PHENCYCLIDINE (PCP 25)</b>			
Phencyclidine	25	4-Hydroxyphencyclidine	12,500
<b>PROPOXYPHENE (PPX 300)</b>			
D-Propoxyphene	300	D-Norpropoxyphene	300
<b>TRICYCLIC ANTIDEPRESSANTS (TCA 1,000)</b>			
Nortriptyline	1,000	Imipramine	400
Nordoxepine	500	Clomipramine	50,000
Trimipramine	3,000	Doxepine	2,000
Amitriptyline	1,500	Maprotiline	2,000
Promazine	3,000	Promethazine	50,000
Desipramine	200	Perphenazine	50,000
Cyclobenzaprine	2,000	Dithiaden	10,000
<b>TRICYCLIC ANTIDEPRESSANTS (TCA 500)</b>			
Nortriptyline	500	Imipramine	200
Nordoxepine	250	Clomipramine	25,000
Trimipramine	1,500	Doxepine	1,000
Amitriptyline	750	Maprotiline	1,000
Promazine	1,500	Promethazine	25,000
Desipramine	100	Perphenazine	25,000
Cyclobenzaprine	1,000	Dithiaden	5,000
<b>TRICYCLIC ANTIDEPRESSANTS (TCA 300)</b>			
Nortriptyline	300	Imipramine	120
Nordoxepine	150	Clomipramine	15,000
Trimipramine	900	Doxepine	600
Amitriptyline	450	Maprotiline	600
Promazine	900	Promethazine	15,000
Desipramine	60	Perphenazine	15,000
Cyclobenzaprine	800	Dithiaden	3,000
<b>TRAMADOL (TML 100)</b>			
n-Desmethyl-cis-tramadol	200	o-Desmethyl-cis-tramadol	10,000
Cis-tramadol	100	Phencyclidine	100,000
Procyclidine	100,000	d,l-O-Desmethyl venlafaxine	50,000
<b>TRAMADOL (TML 200)</b>			
n-Desmethyl-cis-tramadol	400	o-Desmethyl-cis-tramadol	20,000
Cis-tramadol	200	Phencyclidine	200,000
Procyclidine	200,000	d,l-O-Desmethyl venlafaxine	100,000
<b>TRAMADOL (TML 300)</b>			
n-Desmethyl-cis-tramadol	600	o-Desmethyl-cis-tramadol	30,000
Cis-tramadol	300	Phencyclidine	300,000
Procyclidine	300,000	d,l-O-Desmethyl venlafaxine	150,000
<b>TRAMADOL (TML 500)</b>			
n-Desmethyl-cis-tramadol	1,000	o-Desmethyl-cis-tramadol	50,000
Cis-tramadol	500	Phencyclidine	500,000
Procyclidine	500,000	d,l-O-Desmethyl venlafaxine	250,000
<b>KETAMINE (KET 1, 000)</b>			
Ketamine	1,000	Benzphetamine	25,000
Dextromethorphan	2,000	(+) Chlorpheniramine	25,000
Methoxyphenamine	25,000	Clonidine	100,000
d-Norpropoxyphene	25,000	EDDP	50,000
Promazine	25,000	4-Hydroxyphencyclidine	50,000
Promethazine	25,000	Levorphanol	50,000
Pentazocine	25,000	MDE	50,000
Phencyclidine	25,000	Meperidine	25,000
Tetrahydrozoline	500	β-Methamphetamine	50,000
Mephentermine	25,000	l-Methamphetamine	50,000
(1R, 2S) - (-)-Ephedrine	100,000	3,4-Methylenedioxy-methamphetamine (MDMA)	100,000
Disopyramide	25,000	Thionidazine	50,000
<b>KETAMINE (KET 500)</b>			
Ketamine	500	Benzphetamine	12,500

Dextromethorphan	1,000	(+) Chlorpheniramine	12,500
Methoxyphenamine	12,500	Clonidine	50,000
d-Norpropoxyphene	12,500	EDDP	25,000
Promazine	12,500	4-Hydroxyphenacylidine	25,000
Promethazine	12,500	Levorphanol	25,000
Pentazocine	12,500	MDE	25,000
Phencyclidine	12,500	Meperidine	12,500
Tetrahydrozoline	250	d-Methamphetamine	25,000
Mephentermine	12,500	l-Methamphetamine	25,000
(1R, 2S) - (-)-Ephedrine	50,000	3,4-Methylenedioxyamphet- amine (MDMA)	50,000
Disopyramide	12,500	Thioridazine	25,000
<b>KETAMINE (KET 300)</b>			
Ketamine	300	Benzphetamine	6,250
Dextromethorphan	800	(+) Chlorpheniramine	6,250
Methoxyphenamine	8,250	Clonidine	30,000
d-Norpropoxyphene	8,250	EDDP	15,000
Promazine	8,250	4-Hydroxyphenacylidine	15,000
Promethazine	8,250	Levorphanol	15,000
Pentazocine	8,250	MDE	15,000
Phencyclidine	8,250	Meperidine	6,250
Tetrahydrozoline	150	d-Methamphetamine	15,000
Mephentermine	8,250	l-Methamphetamine	15,000
(1R, 2S) - (-)-Ephedrine	30,000	3,4-Methylenedioxyamphet- amine (MDMA)	30,000
Disopyramide	8,250	Thioridazine	15,000
<b>KETAMINE (KET 100)</b>			
Ketamine	100	Benzphetamine	2,000
Dextromethorphan	200	(+) Chlorpheniramine	2,000
Methoxyphenamine	2,000	Clonidine	10,000
d-Norpropoxyphene	2,000	EDDP	5,000
Promazine	2,000	4-Hydroxyphenacylidine	5,000
Promethazine	2,000	Levorphanol	5,000
Pentazocine	2,000	MDE	5,000
Phencyclidine	2,000	Meperidine	2,000
Tetrahydrozoline	50	d-Methamphetamine	5,000
Mephentermine	2,000	l-Methamphetamine	5,000
(1R, 2S) - (-)-Ephedrine	10,000	Thioridazine	5,000
Disopyramide	2,000	3,4-Methylenedioxyamphet- amine (MDMA)	10,000
<b>OXYCODONE (OXY 300)</b>			
Oxycodone	300	Hydromorphone	150,000
Oxymorphone	900	Naloxone	75,000
Levorphanol	15,000	Naltrexone	75,000
Hydrocodone	75,000		
<b>OXYCODONE (OXY 100)</b>			
Oxycodone	100	Hydromorphone	50,000
Oxymorphone	300	Naloxone	25,000
Levorphanol	50,000	Naltrexone	25,000
Hydrocodone	25,000		
<b>COTININE (COT 300)</b>			
(-)-Cotinine	300	(-)-Nicotine	7,500
<b>COTININE (COT 200)</b>			
(-)-Cotinine	200	(-)-Nicotine	5,000
<b>COTININE (COT 100)</b>			
(-)-Cotinine	100	(-)-Nicotine	2,500
<b>COTININE (COT 500)</b>			
(-)-Cotinine	500	(-)-Nicotine	12,500
<b>COTININE (COT 50)</b>			
(-)-Cotinine	50	(-)-Nicotine	1,250
<b>COTININE (COT 10)</b>			
(-)-Cotinine	10	(-)-Nicotine	250
<b>2-ETHYLIDENE-1,5-DIMETHYL-3,3-DIPHENYLPYRROLIDINE (EDDP 300)</b>			

<b>2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)</b>			
300			
<b>2-ETHYLIDENE-1,5-DIMETHYL-3,3-DIPHENYLPYRROLIDINE (EDDP 100)</b>			
100			
<b>2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)</b>			
100			
<b>FENTANYL (FYL 300)</b>			
Alfentanil	>800,000	Buspirone	80,000
Norfentanyl	80	Fentanyl	300
Fenfuramine	150,000	Sufentanyl	150,000
<b>FENTANYL (FYL 200)</b>			
Alfentanil	>800,000	Buspirone	30,000
Fenfuramine	100,000	Fentanyl	200
Norfentanyl	40	Sufentanyl	100,000
<b>FENTANYL (FYL 100)</b>			
Alfentanil	800,000	Buspirone	15,000
Fenfuramine	50,000	Fentanyl	100
Norfentanyl	20	Sufentanyl	50,000
<b>FENTANYL (FYL 20)</b>			
Alfentanil	800,000	Buspirone	15,000
Fenfuramine	50,000	Fentanyl	100
Norfentanyl	20	Sufentanyl	50,000
paliperidone	1,250	Risperidone	5,000
<b>FENTANYL (FYL 10)</b>			
Alfentanil	300,000	Buspirone	8,000
Fenfuramine	25,000	Fentanyl	50
Norfentanyl	10	Sufentanyl	25,000
paliperidone	500	Risperidone	2,500
<b>SYNTHETIC MARIJUANA (K2-50)</b>			
JWH-018 5-Pentanoic acid	50	JWH-073 4-butanoic acid	50
JWH-018 4-Hydroxypentyl	400	JWH-018 5-Hydroxypentyl	500
JWH-073 4-Hydroxybutyl	500		
<b>SYNTHETIC MARIJUANA (K2-30)</b>			
JWH-018 5-Pentanoic acid	30	JWH-073 4-butanoic acid	30
JWH-018 4-Hydroxypentyl	250	JWH-018 5-Hydroxypentyl	300
JWH-073 4-Hydroxybutyl	300		
<b>SYNTHETIC MARIJUANA (K2-25)</b>			
JWH-018 5-Pentanoic acid	25	JWH-073 4-butanoic acid	25
JWH-018 4-Hydroxypentyl	200	JWH-018 5-Hydroxypentyl	250
JWH-073 4-Hydroxybutyl	250		
<b>6-MONOACETYLMORPHINE (6-MAM 10)</b>			
6-Monoacetylmorphine	10	Morphine	100,000
<b>(±) 3, 4-METHYLENEDIOXYAMPHETAMINE (MDA 500)</b>			
(±) 3,4-Methylenedioxy amphetamine	500	Methoxyphenamine	8,000
D,L-Amphetamine sulfate	300	D-Amphetamine	2,000
L-Amphetamine	25,000	Phentermine	1,000
		Maprotiline	50,000
<b>ETHYL-β-D-GLUCURONIDE (ETG 300)</b>			
Ethyl-β-D-Glucuronide	300	Propyl β-D-glucuronide	30,000
Morphine 3β-glucuronide	80,000	Morphine 6β-glucuronide	80,000
Glucuronic Acid	80,000	Ethanol	>100,000
Methanol	>100,000		
<b>ETHYL-β-D-GLUCURONIDE (ETG 500)</b>			
Ethyl-β-D-Glucuronide	500	Propyl β-D-glucuronide	50,000
Morphine 3β-glucuronide	100,000	Morphine 6β-glucuronide	100,000
Glucuronic Acid	100,000	Ethanol	>100,000
Methanol	>100,000		
<b>ETHYL-β-D-GLUCURONIDE (ETG 1,000)</b>			
Ethyl-β-D-Glucuronide	1,000	Propyl β-D-glucuronide	100,000
Morphine 3β-glucuronide	>100,000	Morphine 6β-glucuronide	>100,000
Glucuronic Acid	>100,000	Ethanol	>100,000
Methanol	>100,000		
<b>ETHYL-β-D-GLUCURONIDE (ETG 1,500)</b>			
Ethyl-β-D-Glucuronide	1,500	Propyl β-D-glucuronide	150,000
Morphine 3β-glucuronide	>100,000	Morphine 6β-glucuronide	>100,000

Glucuronic Acid	>100,000	Ethanol	>100,000
Methanol	>100,000		
<b>CLONAZEPAM (CLO 400)</b>			
Clonazepam	400	Flunitrazepam	300
Alprazolam	200	(±) Lorazepam	1,250
a-hydroxyalprazolam	2,000	RS-Lorazepamglucuronide	250
Bromazepam	1,000	Midazolam	5,000
Chlordiazepoxide	1,000	Nitrazepam	200
Clobazam	250	Norchlordiazepoxide	200
Clorazepatedipotassium	800	Nordiazepam	1,000
Delorazepam	1,000	Oxazepam	350
Desalkylfurazepam	250	Temazepam	150
Diazepam	300	Triazolam	5,000
Estazolam	1,250		
<b>CLONAZEPAM (CLO 150)</b>			
Clonazepam	150	Flunitrazepam	120
Alprazolam	75	(±) Lorazepam	500
a-hydroxyalprazolam	750	RS-Lorazepamglucuronide	100
Bromazepam	400	Midazolam	2,000
Chlordiazepoxide	400	Nitrazepam	75
Clobazam	100	Norchlordiazepoxide	75
Clorazepatedipotassium	250	Nordiazepam	400
Delorazepam	400	Oxazepam	130
Desalkylfurazepam	100	Temazepam	60
Diazepam	120	Triazolam	2,000
Estazolam	500		
<b>LYSERGIC ACID DIETHYLAMIDE (LSD 10)</b>			
Lysergic Acid Diethylamide	10		
<b>LYSERGIC ACID DIETHYLAMIDE (LSD 20)</b>			
Lysergic Acid Diethylamide	20		
<b>LYSERGIC ACID DIETHYLAMIDE (LSD 50)</b>			
Lysergic Acid Diethylamide	50		
<b>METHYLPHENIDATE (MPD 300)</b>			
Methylphenidate (Ritalin)	300	Ritalinic Acid	1,000
<b>METHYLPHENIDATE (MPD 150)</b>			
Methylphenidate (Ritalin)	150	Ritalinic Acid	500
<b>METHYLPHENIDATE (MPD 1,000)</b>			
Methylphenidate (Ritalin)	850	Ritalinic Acid	1,000
<b>ZOLPIDEM (ZOL 50)</b>			
Zolpidem	50		
<b>MEPHEDRONE (MEP 500)</b>			
Mephedrone HCl	500	R(+)-Methcathinone HCl	7,500
S(-)-Methcathinone HCl	2,500	3-Fluoromethcathinone HCl	7,500
4-Fluoromethcathinone HCl	1,500	Methoxyphenamine	100,000
<b>MEPHEDRONE (MEP 100)</b>			
Mephedrone HCl	100	R(+)-Methcathinone HCl	1,500
S(-)-Methcathinone HCl	500	3-Fluoromethcathinone HCl	1,500
4-Fluoromethcathinone HCl	300	Methoxyphenamine	100,000
<b>3, 4-METHYLENEDIOXYPYROVALERONE (MDPV 1,000)</b>			
3, 4- methylenedioxy- pyrovalerone	1,000		
<b>3, 4-METHYLENEDIOXYPYROVALERONE (MDPV 500)</b>			
3, 4- methylenedioxy- pyrovalerone	500		
<b>3, 4-METHYLENEDIOXYPYROVALERONE (MDPV 300)</b>			
3, 4- methylenedioxy- pyrovalerone	300		
<b>DIAZEPAM (DIA 300)</b>			
Diazepam	300	Midazolam	6,000
Clobazam	200	Nitrazepam	200
Clonazepam	500	Norchlordiazepoxide	100
Clorazepate dipotassium	500	Nordiazepam	900
Alprazolam	100	Flunitrazepam	200

a-hydroxyalprazolam	1,500	(±) Lorazepam	3,000
Bromazepam	900	RS-Lorazepam glucuronide	200
Chlordiazepoxide	900	Triazolam	3,000
Estazolam	8,000	Temazepam	100
Delorazepam	900	Oxazepam	300
Desalkylflurazepam	200		
<b>DIAZEPAM (DIA 200)</b>			
Diazepam	200	Midazolam	4,000
Clobazam	120	Nitrazepam	120
Clonazepam	300	Norchlordiazepoxide	70
Clorazepate dipotassium	300	Nordiazepam	800
Alprazolam	70	Flunitrazepam	120
a-hydroxyalprazolam	1,000	(±) Lorazepam	2,000
Bromazepam	800	RS-Lorazepam glucuronide	120
Chlordiazepoxide	900	Triazolam	2,000
Estazolam	4,000	Temazepam	70
Delorazepam	900	Oxazepam	200
Desalkylflurazepam	120		
<b>ZOPICLONE (ZOP 300)</b>			
Zopiclone-x-oxide	300	Zopiclone	300
<b>ZOPICLONE (ZOP 50)</b>			
Zopiclone-x-oxide	50	Zopiclone	50
<b>METHCATHINONE (MCAT 500)</b>			
S(-)-Methcathinone HCl	500	R(+)-Methcathinone HCl	1,500
Methoxyphenamine	100,000	3-Fluoromethcathinone HCl	1,500
<b>7-AMINOCLOAZEPAM (7-ACL 300)</b>			
a-hydroxyalprazolam	8,000	Flunitrazepam	3,000
Bromazepam	8,000	RS-Lorazepam glucuronide	2,700
Chlordiazepoxide	8,000	Norchlordiazepoxide	4,500
Clobazam	8,000	Nordiazepam	15,000
Clonazepam	2,400	Temazepam	9,000
Delorazepam	8,000	7-Aminocloazepam	300
Desalkylflurazepam	8,000		
<b>7-AMINOCLOAZEPAM (7-ACL 200)</b>			
a-hydroxyalprazolam	4,000	Flunitrazepam	2,000
Bromazepam	4,000	RS-Lorazepam glucuronide	1,800
Chlordiazepoxide	4,000	Norchlordiazepoxide	3,000
Clobazam	8,000	Nordiazepam	10,000
Clonazepam	1,800	Temazepam	8,000
Delorazepam	4,000	7-Aminocloazepam	200
Desalkylflurazepam	4,000		
<b>7-AMINOCLOAZEPAM (7-ACL 100)</b>			
a-hydroxyalprazolam	2,000	Flunitrazepam	1,000
Bromazepam	2,000	RS-Lorazepam glucuronide	900
Chlordiazepoxide	2,000	Norchlordiazepoxide	1,500
Clobazam	3,000	Nordiazepam	5,000
Clonazepam	800	Temazepam	3,000
Delorazepam	2,000	7-Aminocloazepam	100
Desalkylflurazepam	2,000		
<b>CARFENTANYL (CFYL 500)</b>			
Carfentanyl	500	Fentanyl	100
Sufentanil	50,000	Ramifentanil	10,000
(±)cis-3-Menthylfentanyl	20,000	Butyl fentanyl	150
<b>CARFENTANYL (CFYL 250)</b>			
Carfentanyl	250	Fentanyl	50
Sufentanil	25,000	Ramifentanil	5,000
(±)cis-3-Menthylfentanyl	10,000	Butyl fentanyl	75
<b>CAFFEINE (CAF 1,000)</b>			
Caffeine	1,000		
<b>CATHINE (CAT 150)</b>			
(+)-Norpseudoephedrine HCl (Cathine)	150	(+)-3,4-Methylenedioxyamphetamine (MDA)	100

d/l-Amphetamine	100	p-Hydroxyamphetamine	100
Tryptamine	12,500	Methoxyphenamine	12,500
<b>TROPICAMIDE (TRO 350)</b>			
Tropicamide	350		
<b>ALPRAZOLAM (ALP 100)</b>			
Benzodiazepines	300	Flunitrazepam	200
a-hydroxyalprazolam	1,500	(±) Lorazepam	3,000
Bromazepam	900	RS-Lorazepamglucuronide	200
Chlordiazepoxide	900	Midazolam	8,000
Clobazam	200	Nitrazepam	200
Clonazepam	500	Norchlordiazepoxide	100
Clorazepatedipotassium	500	Nordiazepam	900
Delorazepam	900	Oxazepam	300
Desalkylflurazepam	200	Temazepam	100
Diazepam	300	Triazolam	3,000
Estazolam	8,000	Alprazolam	100
<b>PREGABALIN (PGB 50,000)</b>			
Pregabalin	50,000		
<b>PREGABALIN (PGB 500)</b>			
Pregabalin	500		
<b>ZALEPLON (ZAL 100)</b>			
Zaleplon	100		
<b>CANNABINOL(CNB 500)</b>			
cannabinol	500	Δ <sup>9</sup> -THC	10,000
11-nor-Δ <sup>9</sup> -THC-9 COOH	300		
<b>GABAPENTIN (GAB 2,000)</b>			
Gabapentin	2,000		
<b>TRAZODONE (TZD 200)</b>			
Trazodone	200		
<b>CARISOPRODOL (CAR 2,000)</b>			
Carisoprodol	2,000		
<b>CARISOPRODOL (CAR 1,000)</b>			
Carisoprodol	1,000		
<b>CARISOPRODOL (CAR 500)</b>			
Carisoprodol	500		
<b>AB-PINACA (ABP 10)</b>			
AB-PINACA	10	AB-PINACA 5-Pentanoic	10
AB-PINACA 5-hydroxypentyl	10	AB-FUBINACA	10
AB-PINACA 4-hydroxypentyl	10,000	UR-144 5-Pentanoic	5,000
UR-144 5-hydroxypentyl	10,000	UR-144 4-hydroxypentyl	10,000
APINACA 5-hydroxypentyl	10,000	ADB-PINACA Pentanoic Acid	10
ADB-PINACA N-(5-hydroxypentyl)	30	5-fluoro AB-PINACA N-(4-hydroxypentyl)	30
5-fluoro AB-PINACA	25		
<b>UR-144 (25)</b>			
UR-144 5-Pentanoic acid	25	UR-144 4-hydroxypentyl	10,000
UR-144 5-hydroxypentyl	5000	XLR-11 4-hydroxypentyl	2,000
5-fluoro AB-Pinaca N-(4-hydroxypentyl)	10,000	ADB-PINAC N-(4-hydroxypentyl)	>10,000
AB-PINACA 4-hydroxypentyl	>10,000		
<b>QUETIAPINE (QTP 1,000)</b>			
Quetiapine	1,000	Norquetiapine	10,000
<b>FLUOXETINE (FLX 500)</b>			
Fluoxetine	500		
<b>KRATOM (KRA 300)</b>			
Mitragynine	300	7-hydroxymitragynine	>50,000
<b>TILIDINE (TLD 50)</b>			
Nortilidine	50	Tilidine	100
<b>ALPHA-PYRROLIDINOVALEROPHENONE (α-PVP 2,000)</b>			
Alpha-Pyrrolidinovalerophenone	2,000		
<b>ALPHA-PYRROLIDINOVALEROPHENONE (α-PVP 1,000)</b>			
Alpha-Pyrrolidinovalerophenone	1,000		
<b>ALPHA-PYRROLIDINOVALEROPHENONE (α-PVP 500)</b>			

Alpha-Pyrrolidinovalerophenone	500		
<b>ALPHA-PYRROLIDINOVALEROPHENONE (α-PVP 300)</b>			
Alpha-Pyrrolidinovalerophenone	300		
<b>MESCALINE (MES 100)</b>			
Mescaline	100		
<b>MESCALINE (MES 300)</b>			
Mescaline	300		
<b>PAPAVERINE (PAP 500)</b>			
Papaverine	500	Difunisal	1,000,000
Methortrexate	85,000	Methedrone	500,000
Pragablin	500,000	Phenelzine	8,000
Quinine	4,000		
<b>TAPENTADOL (TAP 1,000)</b>			
3-((1R,2R)-3-(dimethylamino)-1-ethyl-2-methylpropyl)phenol	1,000		
<b>CITALOPRAM (CIT 500)</b>			
Desmethylcitalopram	500		
<b>F-KETAMINE (FKET 1,000)</b>			
2-(2-fluorophenyl)-2-methylamino-cyclohexanone	1,000		
<b>RISPERIDONE (RPD 150)</b>			
Risperidone	150		
<b>SCOPOLAMINE (SCOP 500)</b>			
Scopolamine	500	atropine	3,000
<b>N, N-DIMETHYLTRYPTAMINE (NND 1,000)</b>			
N, N-Dimethyltryptamine	1,000		
<b>MIRTAZAPINE (MTZ 500)</b>			
N-Desmethylmirtazapine	500	Mirtazapine	500
<b>OLANZAPINE (OZP 1,000)</b>			
Olanzapine	1,000		
<b>HYDROMORPHONE (HMO 500)</b>			
Hydromorphone	500	Morphine	200
Codeine	120	Ethylmorphine	120
Hydrocodone	500	Morphine 3-β-D-Glucuronide	250
Levorphanol	2,000	Oxycodone	125,000
Normorphine	125,000	Norcodeine	312,000
Oxymorphone	125,000	Nalorphine	50,000
Thebaine	10,000	Diacetylmorphine (Heroin)	250
8-Monoacetylmorphine	120		
<b>HYDROMORPHONE (HMO 300)</b>			
Hydromorphone	300	Morphine	120
Codeine	75	Ethylmorphine	75
Hydrocodone	300	Morphine 3-β-D-Glucuronide	150
Levorphanol	1,200	Oxycodone	75,000
Normorphine	75,000	Norcodeine	18,700
Oxymorphone	75,000	Nalorphine	30,000
Thebaine	8,000	Diacetylmorphine (Heroin)	150
8-Monoacetylmorphine	75		
<b>HYDROMORPHONE (HMO 250)</b>			
Hydromorphone	250	Morphine	100
Codeine	80	Ethylmorphine	80
Hydrocodone	250	Morphine 3-β-D-Glucuronide	125
Levorphanol	1,000	Oxycodone	62,500
Normorphine	82,500	Norcodeine	15,600
Oxymorphone	82,500	Nalorphine	25,000
Thebaine	5,000	Diacetylmorphine (Heroin)	125
8-Monoacetylmorphine	80		

#### Vliv specifické hmotnosti moči

Patnáct (15) vzorků moči s normálním, vysokým a nízkým rozsahem specifické hmotnosti (1,005-1,045) bylo obohaceno drogami. Koncentrace byla upravena na 50 % pod a 50 % nad cut-off testu. **Multi-Drug Rapid Test** byl testován v duplikátech za použití patnácti vzorků moči bez drog a obohacené moči. Výsledky ukazují, že různé rozsahy specifické hmotnosti moči neovlivňují výsledky testu.

#### Vliv pH moči

pH alikvotního negativního poolu moči bylo upraveno na pH v rozmezí 5 až 9 v přírůstcích po 1 jednotce pH a doplněno drogami. Koncentrace byla upravena na 50 % pod a 50 % nad cut-off testu. Obohacená moč s upraveným pH byla testována pomocí **Multi-Drug Rapid Testu**. Výsledky ukazují, že různé rozsahy pH neovlivňují provedení testu.

#### Zkřížená reaktivita

Byla provedena studie za účelem stanovení zkřížené reaktivity testu se sloučeninami buď v moči bez drog, nebo v moči pozitivní na drogy obsahující výše uvedené související kalibrační látky. Následující sloučeniny nevykazují žádnou zkříženou reaktivitu při testování pomocí Multi-Drug Rapid Test v koncentraci 100µ/ml.

#### Látky nezpůsobující zkříženou reakci

Acetophenetidin	Cortisone	Zomepirac	Quinidine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinine
Acetylsalicylic acid	Deoxycorticosterone	Labetalol	Salicylic acid
Aminopyrine	Dextromethorphan	Loperamide	Serotonin
Amoxicillin	Diclofenac	Meprobamate	Sulfamethazolin
Ampicillin	Diffunisal	Isoxsuprine	Sulindac
I-Ascorbic acid	Digoxin	d,l-Propranolol	Tetracycline
Apomorphine	Diphenhydramine	Nalidixic acid	Tetrahydrocortisone,
Aspartame	Ethyl-p-aminobenzoate	Naproxen	3-acetate
Atropine	β-Estradiol	Niacinamide	Tetrahydrocortisone
Benzilic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrozoline
Benzoic acid	Erythromycin	Norethindrone	Thiamine
Bilirubin	Fenoprofen	Noscapine	Thioridazine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	d,l-Tyrosine
Cannabidiol	Gentisic acid	Oxalic acid	Tolbutamide
Chloral hydrate	Hemoglobin	Oxolinic acid	Triamterene
Chloramphenicol	Hydralazine	Oxymetazoline	Trifluoperazine
Chlorothiazide	Hydrochlorothiazide	Papaverine	Trimethoprim
d,l-Chlorpheniramine	Hydrocortisone	Penicillin-G	d,l-Tryptophan
Chlorpromazine	o-Hydroxyhippuric acid	Perphenazine	Uric acid
Cholesterol	3-Hydroxytyramine	Phenelzine	Verapamil
Clonidine	d,l-Isoproterenol	Prednisone	

#### 【VLASTNOSTI ALKOHOLU】

Detekční limit rychlého testu na alkohol v moči je od 0,02 % do 0,30 % pro přibližnou relativní hladinu alkoholu v krvi. Hraniční hladina alkoholu v moči se může lišit v závislosti na místních předpisech a zákonech. Výsledky testů lze porovnat s referenčními úrovněmi pomocí barevné tabulky na obalu fólie.

#### 【SPECIFICITA TESTU NA ALKOHOL】

Rychlý test na alkohol v moči bude reagovat s methyl, ethyl a allyl alkoholy.

#### 【LÁTKY INTERFERUJÍCÍ S ALKOHOLEM】

Následující látky mohou interferovat s rychlým testem na alkohol v moči při použití jiných vzorků než moči. Jmenované látky se normálně nevyskytují v moči v dostatečném množství, aby interferovaly s testem.

A. Látky, které zlepšují vývoj barvy

• Peroxidázy • Silná oxidační činidla

B. Látky, které inhibují vývoj barvy

• Redukční činidla: kyselina askorbová, kyselina třísllová, pyrogalol, merkaptany a tosyláty, kyselina šťavelová, kyselina močová





• Bilirubin • L-dopa

• L-methylidopa • Metampyrin

#### 【BIBLIOGRAFIE】

1. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735.
2. B. Cody, J.T., "Specimen Adulteration in drug urinalysis. Forensic Sci. Rev., 1990, 2:63.
3. C. Tsai, S.C. et.al., J. Anal. Toxicol. 1998; 22 (6): 474
4. Hawks RL, CN Chiang. *Urine Testing for Drugs of Abuse*. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

#### Seznam symbolů

	Prostudujte návod k použití		Počet testů v soupravě		Autorizovaný zástupce
	Pouze pro in vitro diagnostiku		Spotřebujte do		Nepoužívat opakovaně
	30°C Skladujte mezi 2-30°C		Číslo šarže		Katalog. číslo
	Nepoužívejte v případě poškození obalu		výrobce		



**Hangzhou AllTest Biotech Co.,Ltd.**  
#550, Yinhai Street  
Hangzhou Economic & Technological Development Area  
Hangzhou, 310018 P.R. China  
Web: www.alltests.com.cn Email: info@alltests.com.cn



**MedNet GmbH**  
Borkstrasse 10  
48163 Muenster  
Germany

Number: 145177407  
Effective Date: 2023-03-20

Český překlad :2023-11-07/VE